



Method 5.5 – Syrup and remelt: ICUMSA 420 colour and turbidity

1. Rationale

This method is applicable to factory syrups and remelt and is used to determine the ICUMSA colour and/or turbidity of a solution of the sample at pH 7.00.

2. Principle

A diluted sample is divided in two portions. One portion is filtered through a membrane filter (the type of which is specified by ICUMSA) to remove turbidity. The pH of both portions is adjusted to 7.00 ± 0.02 using basic or acidic solutions. The Brix and absorbance of the solutions are measured at a wavelength of 420 nm and the ICUMSA values of the solutions calculated. The colour of the sample is the ICUMSA value of the filtered portion and the turbidity is the difference between the ICUMSA values of the filtered and unfiltered portions.

3. Definitions

3.1 Transmittance of a solution

If I_1 represents the radiant energy incident upon the first surface of the solution, and I_2 represents the radiant energy leaving the second surface of the solution. Then:

$$T = \frac{I_2}{I_1} = \text{transmittance of the solution}$$

and $100 \times T = \text{percentage transmittance}$

3.2 Transmittancy

Let T_{soln} represent the transmittance of a cell containing the solution and let T_{solv} represent the transmittance of the same cell containing the pure solvent. Then:

$$T_s = \frac{T_{\text{soln}}}{T_{\text{solv}}} = \text{transmittancy of the solution}$$

3.3 Absorbancy (extinction) measured in absorbance units (AU)

$$A_s = -\log_{10} T_s = \text{absorbancy of the solution}$$

3.4 Absorbancy index (extinction index)

Let b represent the length (mm) of the absorbing path between the boundary layers of the solution and let c represent the concentration (g/cm^3) of the sugar solution. Then:

$$A_i = \frac{A_s}{bc} = \text{absorbancy index of the solution.}$$

3.5 ICUMSA Colour

The value of the absorbancy index multiplied by 10 000 is reported as the ICUMSA Colour of the solution and the resulting value is expressed in ICUMSA Units (IU). Since the wavelength at which the determination of colour in solution is used is set at 420 nm the value is designated as being the ICUMSA 420 Colour.

4. Apparatus

- 4.1 **Spectrophotometer** capable of light transmission measurements at a wavelength of 420 nm with the narrowest practical bandwidth, *e.g.* ± 10 nm
- 4.2 **Optical glass cell**: 10 mm
- 4.3 **Membrane filters**: cellulose nitrate filters, 0.45 μm pore size, 50 or 47 mm ϕ
- 4.4 **Buchner funnel** or magnetic vacuum filtration funnel: 50-65 mm ϕ
- 4.5 **Buchner flask** and rubber bung: 500 cm^3
- 4.6 **pH meter** capable of measuring to 0.01 units
- 4.7 **Refractometer**
- 4.8 **Magnetic stirrer** with stirrer bar
- 4.9 **Laboratory balance** readable to 0.01 g
- 4.10 **Beakers**: 100, 250 and 1 000 cm^3
- 4.11 **Measuring cylinder**: 100 cm^3
- 4.12 **Volumetric flasks**: 2 \times 200, 2 \times 1 000 cm^3
- 4.13 **Pipettes**: 2 \times 10 cm^3

5. Reagents

5.1 Hydrochloric acid solution (1 M)

Hydrochloric acid (HCl, 32%) is a corrosive acid and contact with the skin, eyes and through inhalation must be avoided. Work in a fume cupboard while wearing gloves and safety glasses.

Measure 98 cm^3 concentrated hydrochloric acid and add to approximately 700 cm^3 distilled water in a 1 000 cm^3 beaker. Always add the acid to the water and not the other way around. The acid dilution is exothermic and the solution will therefore heat. Allow the solution to cool down, transfer to a 1 000 cm^3 volumetric flask and make to the mark.

5.2 Hydrochloric acid solution (0.05 M)

Pipette 10 cm^3 of the 1 M hydrochloric acid solution into a 200 cm^3 volumetric flask and make to the mark with distilled water.

5.3 Sodium hydroxide solution (1 M)

Sodium hydroxide (NaOH) is a corrosive base and contact with the skin and eyes should be avoided. Wear gloves and safety glasses during use.

Weigh 40.0 g sodium hydroxide pellets and dissolve in some distilled water. This dissolution is exothermic and the solution will therefore heat. Allow the solution to cool and dilute to 1 000 cm^3 in a volumetric flask.

5.4 Sodium hydroxide solution (0.05 M)

Pipette 10 cm³ of the 1 M sodium hydroxide solution into a 200 cm³ volumetric flask and make to the mark with distilled water.

6. Procedure

6.1 Calibration of pH meter

Following the manufacturer's directions, calibrate the pH meter using the 4.00 and 7.00 pH buffer solutions (compensated for a temperature different from 20°C) while stirring at a constant rate. Calibrations should be done at the beginning of each day or shift using fresh buffer solutions only. The buffer solutions should be at room temperature.

6.2 Sample analysis

Allow the clear juice to cool to room temperature. Mix the sample thoroughly. Determine the approximate Brix of the unfiltered solution and calculate the mass of juice needed to prepare a 5°Bx solution using the formula below.

$$\text{Mass juice (g)} = \frac{\text{total mass (g)}}{\text{original Brix (° Bx)}} \times \text{final Brix (° Bx)}$$

Dilute this mass of well-mixed sample to 100.00 g using distilled water.

For analysis of colour and turbidity, divide the solution into two parts and filter the one through a 0.45 µm cellulose nitrate membrane under vacuum into a clean, dry Buchner flask. For the analysis of colour only filter the whole solution through a 0.45 µm membrane under vacuum into a clean, dry Buchner flask.

Transfer the solutions into 50 or 100 cm³ beakers. Stir the solutions on the magnetic stirrer and adjust the pH to 7.00 ± 0.02 using either hydrochloric acid (0.05 M) to bring the pH down or sodium hydroxide (0.05 M) to bring the pH up. Allow ample time for the pH reading to stabilize (1 minute).

Measure the absorbancies of the solutions in a 10 mm cell using the spectrophotometer at 420 nm against distilled water as a reference. Use a 5 mm cell if the absorbance readings are unstable. Also measure the Brix of the solutions in a refractometer.

7. Calculations

Use the Brix readings to obtain the concentration of total solids in g/cm³ using the formula indicated below. Determine the ICUMSA colour of the filtered and unfiltered solutions.

$$\text{Total solids (g/cm}^3\text{)} = \frac{\text{Brix} \times (\text{a} + \text{b} \times \text{Brix} + \text{c} \times \text{Brix}^2 + \text{d} \times \text{Brix}^3 + \text{e} \times \text{Brix}^4)}{100}$$

where a	=	0.9971843
b	=	3.85738 × 10 ⁻³
c	=	1.254916 × 10 ⁻⁵
d	=	8.125659 × 10 ⁻⁸
e	=	5.611455 × 10 ⁻¹⁰

$$\text{ICUMSA 420 Colour} = \frac{A_s \times 10\,000}{bc}$$

where A_s \equiv absorbance at 420 nm (AU)
 b \equiv cell length (mm)
 c \equiv concentration of total solids (g/cm³)

Report results in ICUMSA Units (IU) to the nearest 10 units.

8. Example

Obtain the concentration of total solids in g/cm³ from the Brix reading using the equation in 7.

Unfiltered solution:

Brix = 5.0°Bx
 Concentration of total solids = 0.051 g/cm³
 Absorbance at 420 nm = 0.985 AU

$$\begin{aligned} \text{ICUMSA 420 colour} &= \frac{0.985 \text{ AU} \times 10\,000}{10 \text{ mm} \times 0.051 \text{ g/cm}^3} \\ &= 19\,314 \text{ IU} \end{aligned}$$

Filtered solution:

Brix = 5.0°Bx
 Concentration of total solids = 0.051 g/cm³
 Absorbance at 420 nm = 1.042 AU

$$\begin{aligned} \text{ICUMSA 420 colour} &= \frac{1.042 \text{ AU} \times 10\,000}{10 \text{ mm} \times 0.051 \text{ g/cm}^3} \\ &= 20\,431 \text{ IU} \end{aligned}$$

$$\begin{aligned} \text{Turbidity} &= \text{Unfiltered Colour} - \text{Filtered Colour} \\ &= (20\,431 - 19\,314) \text{ IU} \\ &= 1\,117 \text{ IU} \end{aligned}$$

Report as 1 117 IU

9. References

SASTA (1985). *Laboratory Manual for South African Sugar Factories*. 3rd Edition, p 287 - 290.

SMRI (1997). Determination of the colour and turbidity of syrup. *SMRI Test Methods*, TM024.